GRIPKIT

BY WEISS ROBOTICS

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DOOSAN

1 Introduction

GRIPKIT CR PLUS is Weiss Robotics' intelligent gripping solution for robot arms of Doosan Robotics. GRIPKIT CR PLUS is available with different gripping modules and is mechanically and electrically fully compatible to the robot arms of Doosan Robotics.

This operating manual describes the installation, initial setup and operation of the GRIPKIT CR PLUS. It does not replace the operating manual of the used gripper which can be found Online at www.weiss-robotics.com/gripkit-cr-plus.

1.1 Related documents

Related documents for operating the gripping module, technical drawings and 3D models are available at www.weiss-robotics.com.

Please refer to our General Terms and Conditions for additional information on the warranty, available at www.weiss-robotics.com/gtc.

1.2 Target audience

This manual is intended for integrators and operators of robot applications and all persons maintaining and working with this product. It is advised to keep this and related documents available to the staff at all times. Please read this document carefully and observe in particular the safety and warning notices in this document.

1.3 Notation and symbols

For a better understanding, the following symbols are used in this manual.

Functional or safety-relevant information. Non-compliance may endanger the safety of personnel and the system, damage the device or impair its function.

Additional information for a better understanding of the described facts.

Reference for further information.

2 Basic safety instructions

2.1 Intended use

GRIPKIT CR PLUS was developed for reliably gripping and holding of work pieces or other objects and is intended for mounting on a robot arm. The requirements of the applicable directives as well as the assembly and operation instructions in this document must be observed and adhered to. GRIPKIT CR PLUS may only be used within its defined operating limits and only for industrial applications.

Any other use is considered to be improper, e.g. if the gripping module is used for pressing, curving, lifting, punching or as a clamping or guidance aid for other tools. The manufacturer is not liable for any damage resulting from this improper use.

2.2 Environmental and operating conditions

GRIPKIT CR PLUS and the included gripping module may only be used in the context of its defined application parameters. Make sure, that the module and the gripper fingers are at sufficient size for the application. The environment has to be clean and the ambient temperature must not exceed the specifications in this document. In addition, the environment must be free from splash water and vapors as well as from abrasion or processing dust, except the used gripping module is especially designed for contaminated environments.

2.3 Controlled production

GRIPKIT CR PLUS represents the state of the art and the recognized safety rules at the time of delivery. However, it can present risks if, for example:

- GRIPKIT CR PLUS is not used in accordance with its intended purpose
- GRIPKIT CR PLUS is modified, or not installed or maintained properly
- The EC Machinery Directive, the VDE directives, the safety and accident-prevention regulations valid at the usage site, or the safety and installation notes are not observed

2.3.1 Protective equipment

Provide protective equipment per EC Machinery Directive where required.

2.3.2 Constructional changes, attachments or modifications

Additional drill holes, threads or attachments that are not offered as accessories by

Weiss Robotics may be attached only with written permission of Weiss Robotics.

2.3.3 Special standards

The following harmonized standards are adhered to:

- RF interference, noise and emission according to IEC/EN 61000-6-3
- Fast transients (bursts) on supply and I/O lines according to IEC/EN 61000-4-4
- RF power input on supply and I/O lines according to IEC/EN 61000-4-6
- RF irradiation according to IEC/EN 61000-4-3
- Interference emission according to IEC/EN 61000-6-4 Class A
- Magnetic field with energy frequency according to IEC/EN 61000-4-8
- Electrostatic discharge according to IEC/EN 61000-4-2
- IO-Link communication standard according to IEC 61131-9

2.4 Personnel qualification

The assembly, initial commissioning, maintenance, and repair of the GRIPKIT CR PLUS may be performed only by trained specialist personnel.

Every person called upon by the operator to work on the GRIPKIT CR PLUS must have read and understood the complete operating manual. This applies particularly to occasional personnel such as maintenance personnel.

2.5 Safety-conscious working

Observe the safety and accident-prevention regulations valid at the usage site.

Do not reach into the open mechanism or the movement area of the gripping module. De-energize and set the gripping module pressure-free before performing any work on it.

3 Warranty

The warranty is 12 month from the ex-works delivery date in accordance with the intended use in one shift operation or 20 million gripping cycles and observing the prescribed maintenance and lubrication intervals. Maintenance and lubrication intervals can be found in the operation manual of the included gripping module. Workpiece contacting parts and wearing parts are not covered by this warranty. Refer to our Terms and Conditions for details. GRIPKIT CR PLUS is considered to be defective, if its basic function gripping is no longer given.


4 Model variants

GRIPKIT CR PLUS is available in different model variants.

GRIPKIT CR PLUS	S	L
Gripping module (Part No.)	CRG 30-050 (5010016)	CRG 200-085 (5010017)
Active mounting plate (Part No.)	SF-ISO50-CRG30 (5020061)	SF-ISO50-CRG200 (5020062)

5 Assembly

The gripping module is mounted to the robot arm using the active mounting plate contained in the GRIPKIT CR PLUS. The necessary mounting material (screws, pins, etc.) is included in the GRIPKIT CR PLUS.

Only work on the electrical system of the robot when it is switched off.


Do not kink the cable at the connector! Damage to the connector possible!

With GRIPKIT CR PLUS S: Install the enclosed spacer plates between the flange and the robot. Use the long screws with a diameter of 6 mm for this.

6 Technical specification

Pin	Function
1	IN2
2	IN1
3	INO
4	Reserved
5	Reserved
6	GND
7	+24 V/ 1.5 A IN
8	OUT2
9	OUT 1
10	OUT 0

Power Consumption

Model	Holding at 25%	Holding at 100%	Peak
GRIPKIT - CR PLUS S	0,15 A	0,24 A	0,51
GRIPKIT - CR PLUS L	0,24 A	0,48 A	1,87

6.1 Input Circuit

The inputs of the gripper module support both NPN and PNP signals. The passive state is detected when the input is open (inactive). The active state is detected when the input is actively pulled to either 0V (NPN) or 24V (PNP).

Input	Function GRIPKIT CR PLUS
INO	GRIP active: grab workpiece passive: workpiece released
IN1	INDEX.0 selects the grip to execute
IN2	INDEX.1 selects the grip to execute

Grip	IN1	IN2
0	p	p
1	a	p
2	a	p
3	a	a

6.2 Output circuit

The gripping state is output via the outputs of the gripping module. The outputs are designed as a push-pull circuit

Out 2	OUT1	OUT0	Gripping state
0 V	0 V	0 V	wait for command (IDLE)
0 V	0 V	24 V	release workpiece (RELEASED)
0 V	24 V	0 V	no workpiece detected (NO PART)
0 V	24 V	24 V	holding workpiece (HOLDING)
24 V	0 oder 24 V	0 oder 24 V	Error

7 Scope of delivery

The following components are included in the GRIPKIT CR PLUS package for UR:

- Gripping module (see chapter „Model variants“)
- Active mounting plate (see chapter „Model variants“)
- Connection kit GRIPKIT CR EASY/PLUS (Part No. 5070018)
- For GRIPKIT CR PLUS S: two spacer plates (Part No. 5020063)
- USB flash drive with software (Part No. 5090008)
- Quickstart GRIPKIT CR PLUS (Part No. 5080030)

Additional Accessories

GRIPKIT CR PLUS does not contain any fingers, since those are usually highly application-dependent.

Order matching fingers from our standard product range separately.

7.1 Content of the included USB flash drive

The included USB flash drive contains the software as well as sample programs and the documentation of the gripping module.

8 Technical data

GRIPKIT CR PLUS	S	L
Gripper type	parallel	
Gripping force	15 - 30 N	75 - 200 N
Workpiece weight (form fit/force fit)	2.7 kg 0.15 kg	4.3 kg 0.15 kg
Full stroke	50 mm	85 mm
Current consumption	350 mA	850 mA
Weight	0.46 kg	1.3 kg
Parameterization	4 grips (position, force)	


Exceeding the specified technical data may reduce the life of the gripping module or permanently damage it. If in doubt, please clarify your application with your responsible sales partner.

9 System requirements


GRIPKIT CR PLUS is compatible with the Doosan Cobots type M0609, M1509, M1013 and M0617. For mounting the GRIPKIT CR PLUS size S, additional spacer plates (Part No. 5020063) are required. The spacer plates are included.

10 Device Configurator


To configure your gripper module, you can alternatively use the "Device Configurator" program. You will find this program on the enclosed USB flash drive. Do the installation on a desktop device and open the program. Connect the Gripper to your computer with the enclosed USB cable.

Start the "Device Configurator" program and connect the gripper via "Connect".

Next, activate the control of the gripper and carry out a reference run. Press "Activate" and then "Reference".


ENGLISH
DOOSAN

Now the gripper is ready for use. You can test the adjusted grips 0-3 by pressing "Grip" and "Release". The respective grip parameters are defined with "Limit value no part" and "Limit value part release".

10.1 Grip direction

The gripping direction is determined by the relation of the position limit values to one another. If the „No Part Limit“ is larger than the „Release Limit“, the gripping module performs an internal grip, if the „Release Limit“ is larger than the „No Part Limit“, an external grip will be performed.

External and internal gripping example

10.1.1 Gripping force

GRIPKIT CR PLUS has an integrated gripping force control. The gripping force can be specified as a percentage of the nominal gripping force of the gripping module. Possible values are 25 to 100 % of the nominal gripping force.

10.2 Define a grip

To define a grip, go to the "Gripping Settings" tab. Select the parameters for "Limit value no part" and "Limit value part release" and confirm the entries with "Apply".

When gripping and releasing, electrical gripping modules stop at the „No Part Limit“ and the „Release Limit“.

You must then specify the gripping force as a percentage of the nominal force. The nominal force of the gripping module used can be found in chapter 8 Technical data. With the parameters "Override gripping" and "Override enable" you can readjust the gripping speed. Both values are set to 100 % by default. When gripping, this value refers to the optimal gripping speed calculated by the gripper depending on the part size and gripping force. You can set values from 1-200 % here. When releasing, 100 % is the absolute maximum speed of the gripper. Values of 1-100 % are possible here.

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
10.3 Set up light ring

The gripping modules have a glare-free LED light ring that can be used to visualize various states of the application. Eight visualization patterns can be parameterized for this purpose. The animation, colour and speed can be changed. To do this, open the "LED" tab and select the various gripping states via "Animation". Confirm the input with "Activate".


10.4 Diagnostics


Under the "Diagnostics" tab you will find the current data on the status of the connected gripping module.

11 Installing

To install the software, first switch on the robot controller. Then the robot must be switched on using the enclosed teach pendant. After starting up the Teach Pendant, open the "Status" button.


Select "Servo On" in the next step to activate the robot and test all functions.


ENGLISH
DOOSAN

To integrate the program, the "Task Writer" located on the menu must be opened. Here you can either continue working on an already existing project or create a new project. Then confirm your selected project with "OK".


Then open the cover of your controller and insert the enclosed USB flash drive into the USB port. To import the program you need the "Custom Code" function. Insert this function between the "global variables" function and the "main sub" function to make the best use of them. Select the "Custom Code" function and switch to the "Property" tab.


ENGLISH
DOOSAN

Then activate "Script file" to be able to use the Open button to the right. This makes all .txt files of the USB flash drive visible.

In the next step select the document "gripkit_plus.txt" and click on "Open".

To confirm this change, click on "Confirm" in the top right-hand corner (the button changes to a grey state when pressed). The gripper can now be used in the program.

Always save your changes directly to prevent data loss.

1. switch to the "Play" tab to automatically save the program.


or

2. press the Settings icon in the upper left corner of the Teach Pendant and save your program using the drop-down menu.

11.1 Testing the gripping functions

Switch to the "Status" button and press the "Servo On" button to activate the robot and test all functions.

Now you can control the gripper via the Flange Digital Outputs. To perform a gripping operation, activate the "Flange Digital Output" number 1.

The gripper status can then be checked via the "Flange Digital Inputs". Switch off "Flange Digital Output" number 1 again to enable it.

11.2 Program features

11.2.1 Grip (Syntax: gripkit_grip(grip_index))

This function ensures that the gripper grips. To execute a grip, the "Grip" function requires a value between 1 and 4 for the grip to be executed. The function also returns a return value. This indicates whether the gripper has gripped a workpiece, has not found a workpiece or an error has occurred. The return values are listed in the following table:

Parameter: grip_index: value between 1 and 4

Return	Outcome
2	No Part
3	Holding
>=4	Error

If the gripper is already in the "Holding" or "No Part" state, a release must first be performed before it can be gripped again

The gripper is initialized via the first grip in the program, before the first grip is the state of the gripper "Idle".

11.2.2 Release (Syntax: gripkit_release(grip_index))

This function opens the gripper. It requires a parameter value between 1 and 4 for the release to be executed. The function also returns a return value that indicates whether the gripper is open or whether an error has occurred. The return values are listed in the following table:

Parameter: grip_index: value between 1 and 4

Return	Outcome
1	Released
>=4	Error

If the gripper is already in the "Released" state, a handle must first be executed before it can be re-released.

11.2.3 Gripstate

This function can be used to query the current gripper condition.

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11.3 Programming

For programming, the USB stick containing the program must first be removed. Then connect a USB keyboard to the robot controller.

The following options are available for inserting commands into the program:

1. use the "Custom Code" function There you can transfer the return values of the functions to variables already existing in the program.
2. call the commands directly via an If function. Please note, however, that it is not possible to transfer the return values to a variable, but only to compare values.

Then insert the required command into your program and select "Properties" tab. The following must be observed when inserting:

- With the "Custom Code" function, click directly into the "Singleline" line, or with several commands into the box below, and enter the desired function e.g. assign the return value to a variable defined in the program: `x =gripkit_gripstate()`
- With the If function you click directly into the input line and enter the desired values using a keyboard: e.g. `gripkit_release(4) == 1`

The screenshot displays the 'Task Writer' application. On the left, a 'Tools' sidebar includes icons for Multi-Select, Copy, Cut, Paste, Delete, Row Up, and Row Down. The central 'Task List' pane contains the following tasks:

- 001 GlobalVariables
- 002 CustomCode (#...
- 003 MainSub (Task Vel. 250.0, Acc. 1.0, Jo...
- 004 CustomCode (gripkit_grip(1))
- 005 Define (x)
- 006 CustomCode (x=gripkit_gripstate())
- 007 If (x == 2 or x == 3)
- 008 CustomCode (gripkit_release(3))
- 009 Else
- 010 End If
- 011 EndMainSub

The right-hand pane is divided into 'Command', 'Property', and 'Play' tabs. The 'Custom Code' tab is active, showing a text input field with the code `x=gripkit_gripstate()` and a 'Confirm' button. Below it, the 'Script file' section has an 'Open' button.

To start the program and save it automatically, switch to the "Play" tab. Switch on the "Real Mode", so that the tool outputs can be set. Then press the play icon.

The gripper is initialized via the first grip in the program. That means, before the first grip is the state of the gripper idle. From this it follows that at a release command before a gripper command the gripper is not released.

12 Uninstall

To uninstall the program, delete the "Custom Code" function used for the installation. You should also delete all "Custom Code" and If functions in which the commands described here will be used.

13 Translation of the original declaration of conformity

In terms of the EC Machinery Directive 2006/42/EG, Annex II, Part B

Manufacturer:

Weiss Robotics GmbH & Co. KG
Karl-Heinrich-Kaeferle-Str. 8
71640 Ludwigsburg, Germany

Distributor:

Weiss Robotics GmbH & Co. KG
Karl-Heinrich-Kaeferle-Str. 8
71640 Ludwigsburg, Germany

We hereby declare that the following products

GRIPKIT CR PLUS S, GRIPKIT CR PLUS L

meet the applicable basic requirements of the Machinery Directive (2006/42/EC). The incomplete machine may not be put into operation until conformity of the machine into which the incomplete machine is to be installed with the provisions of the Machinery Directive (2006/42/EC) is confirmed.


Applied harmonized standards, especially:

- EN ISO 12100-1 Safety of machines – Basic concepts, general principles for design – Part 1: Basic terminology, methodology
- EN ISO 12100-2 Safety of machines – Basic concepts, general principles for design – Part 2: Technical principles

The manufacturer agrees to forward the special technical documents for the incomplete machine to state offices on demand. The special technical documents according to Annex VII, Part B, belonging to the incomplete machine have been created.

Person responsible for documentation: Dr.-Ing. K. Weiß, Tel.: +49(0)7141/94702-0

Place, Date, Signature: Ludwigsburg, 7. February 2020

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The specification in this document are subject to change without prior notice for the purpose of product improvement. All trademarks are the property of their respective owners

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Details of the signer: Weiss Robotics GmbH & Co. KG, Management